

E.S.I.T Compatibility Chart



								
Scan Body	Analog	Final Screw	Ti Abutment	Hybrid Zirconia Ti Abutment	Full contour Zirconia Ti Abutment	PFM	3D model	3D model
ESIT Implant			Compatible Implant Company			Size		
<i>ESIT 1</i>			<i>3I Biomet Certain ®*</i>			<i>3.4</i>		
						<i>4.1</i>		
						<i>5.0</i>		
<i>ESIT 2</i>			<i>Noble Active ®*</i>			<i>3.5</i>		
						<i>4.3</i>		
<i>ESIT 3</i>			<i>Noble Replace Select ®*</i>			<i>3.5</i>		
						<i>4.3</i>		
						<i>5</i>		
<i>ESIT 4</i>			<i>Straumann Bonelevel ®*</i>			<i>3.3</i>		
						<i>4.1</i>		
<i>ESIT 5</i>			<i>Osstem ®*</i>			<i>Mini</i>		
						<i>Regular</i>		
<i>ESIT 6</i>			<i>Dentium ®*</i>			<i>4.5</i>		
						<i>4.8</i>		
						<i>6.5</i>		
<i>ESIT 6</i>			<i>Straumann Tissue level ®*</i>			<i>4.8</i>		
						<i>6.5</i>		
<i>ESIT 7</i>			<i>Zimmer ®**</i>			<i>3.5</i>		
						<i>4.5</i>		
						<i>5.7</i>		
<i>ESIT 7</i>			<i>Implant Direct Legacy 3 ®*</i>			<i>3.7 / 4.2</i>		
						<i>4.7 / 5.2</i>		
						<i>5.7</i>		
<i>ESIT 8</i>			<i>I.T.I ®*</i>			<i>3.5</i>		
						<i>4.8</i>		
						<i>6.5</i>		
<i>ESIT 9</i>			<i>Frialit Xive ®*</i>			<i>3.4</i>		
						<i>3.8</i>		
<i>ESIT 10</i>			<i>Dio ®*</i>			<i>3.8</i>		
						<i>4.5</i>		
<i>ESIT 11</i>			<i>Astra Tech bone level ®*</i>			<i>3.0</i>		
						<i>3.5 / 4.0</i>		
						<i>4.5 / 5.0</i>		
<i>ESIT 13</i>			<i>MIS ®*</i>			<i>3.7</i>		
						<i>4.2</i>		
						<i>5.0</i>		

*Compatible companies are not trademarks of Barth Dental Laboratories